



The Abdus Salam
**International Centre
for Theoretical Physics**



2495-09

Joint ICTP-IAEA Workshop on Nuclear Data for Analytical Applications

21 - 25 October 2013

Exercises on PIGE for the Thursday session

Massimo Chiari
*I.N.F.N
Florence*

Exercises on PIGE for the Thursday session

Exercise n. 1

Perform the energy calibration of a gamma-ray detector using the spectrum collected with a ^{60}Co radioactive source (open spectrum Ecalib_60.asc in Spectra_A\ directory). Save the calibration.

Exercise n. 2

Determine the full peak absolute efficiency of the gamma-ray detector using the spectrum collected with a ^{152}Eu radioactive source (open spectrum Efficiency_152Eu.asc in Spectra_A\ directory).

- i) Looking at the “Gamma emission” table in the Eu-152_tables.pdf file select the more intense gamma-rays (i.e. intensity > 1 %) and calculate the peak areas for those lines. If necessary refine the energy calibration.
- ii) Find the ^{152}Eu life-time. Calculate the present source activity knowing that the initial activity was 370 kBq (uncertainty 5%), the time elapsed till the measurement is 10 years and the measurement lasted 4000 s.
- iii) With a spreadsheet calculate and plot the detector full peak absolute efficiency as a function of gamma-ray energy. Calculate the uncertainty.
- iv) Fit the measured efficiency with a polynomial function in the argument of $1/E_{\text{gamma}}$

Exercise 3

Perform the energy calibration of a gamma-ray detector using the PIGE spectrum collected by bombarding a thick Al target with proton (open spectrum Al.asc in Spectra_B\ directory). Save the calibration.

Exercise 4

Calculate the Na (actually Na_2O) content in archeological glass by PIGE measurements with 3 MeV protons in an external beam set-up (this is an exercise on the application of PIGE thick target elemental analysis with comparison to standards).

- i) Using certified thick standards of glass calculate the relationship between the Na_2O percent by weight and the yield of the 441 keV gamma-ray line from the $^{23}\text{Na}(p,p'\gamma)^{23}\text{Na}$ reaction (open spectra multicomp_glass_1412.asc, soda_lime_621.asc and trace_el_glass_610.asc in Spectra_B\ directory). Use the

values reported under the “Norm. factor” column in the Guide to spectra.xls file to normalize the 441 keV peak areas and calculate the yields (and the uncertainties). As a zero-order approach, neglect the correction for the stopping powers. Plot the data and fit them with a line.

- ii) Check the possible contribution from beam induced gamma-ray background to the 441 keV gamma-ray line and eventually calculate it (open spectrum background_beam.asc in Spectra_B\ directory).
- iii) Calculate the yield (and the uncertainty) of the 441 keV gamma-ray line from the $^{23}\text{Na}(p,p'\text{g})^{23}\text{Na}$ reaction from the PIGE spectra of the 11 unknown samples (open spectra red*.asc, green*.asc, blue*.asc, azure*.asc and yellow*.asc in Spectra_B\ directory).
- iv) Calculate the Na₂O percent by weight (and the uncertainty) in the unknown samples using the yields from step iii) and the measurement on the standards from step i).

Exercise 5

Implement the correction for the stopping power in the different matrices in the analysis of PIGE measurements of glass standards.

- i) Using SRIM calculate the actual proton beam energy on the target knowing that the proton beam energy in vacuum is 3 MeV and the beam passes through a 7.5 micron Kapton extraction window and a path of 10 mm of He before striking the target.
- ii) Determine the “E1/2” energy using the thick target yield of the $^{23}\text{Na}(p,p'\text{g})^{23}\text{Na}$ reaction ($E_{\text{gamma}} = 441 \text{ keV}$) retrieved from IBANDL.
- iii) Calculate the stopping powers (for example in $\text{keV}/\mu\text{g}/\text{cm}^2$) for protons of “E1/2” energy in the different matrices of the glass standards, using the composition shown in the Glass standards.xls file in Spectra_B\ directory and converting percent in weight of the oxides to atomic percent of the elements (if needed use the weighttoatomicconversion.xls calculator).
- iv) Replot the relationship between the Na₂O percent by weight of the glass standards and the yield of the 441 keV gamma-ray line from the $^{23}\text{Na}(p,p'\text{g})^{23}\text{Na}$ reaction upon correction for the stopping power.

Exercise 6

Using the composition of the multicomponent glass standard (NIST SRM 1412) and the

table of reactions with prompt emission of gamma-rays induced by protons, identify the gamma-ray lines and the emitting isotopes in the PIGE spectrum of multicomponent glass standard (open spectrum multicomp_glass_1412.asc in Spectra_B\ directory).

Exercise 7

Identify the gamma-ray lines and the emitting isotopes in the spectrum of natural background radiation (open spectrum nat_background.asc in Spectra_A\ directory).