




INTERNATIONAL ATOMIC ENERGY AGENCY  
UNITED NATIONS EDUCATIONAL, SCIENTIFIC AND CULTURAL ORGANIZATION



INTERNATIONAL CENTRE FOR THEORETICAL PHYSICS  
34100 TRIESTE (ITALY) - P.O. B. 586 - MIRAMARE - STRADA COSTIERA 11 - TELEPHONES: 224 411/2/3/4/5/6  
CABLE: CENTRATON - TELEX 480392-1



SMR/115 - 23

WINTER COLLEGE ON LASERS, ATOMIC AND MOLECULAR PHYSICS

(21 January - 22 March 1985)

A COMPACT AND RELIABLE Nd:YAG LASER  
FOR OPHTHALMIC MICROSURGERY

F.DOCCHIO  
L.PALLARO  
O.SVELTO

Centro per l'Elettronica Quantistica  
e Strumentazione Elettronica del CNR  
Istituto di Fisica-Politecnico  
Milan, Italy

M. GUERINI  
Elettronica Valseriana, Laser Div.  
V. Bombardieri, 27  
Fiorano al S., Italy

A COMPACT AND RELIABLE Nd:YAG LASER FOR OPHTHALMIC MICROSURGERY

F. DOCCHIO, L. PALLARO & O. SVELTO  
Centro per l'Elettronica Quantistica  
e Strumentazione Elettronica del CNR  
Istituto di Fisica del Politecnico  
Piazza Leonardo da Vinci, 32  
20133 Milano  
Italy

These are preliminary lecture notes, intended only for distribution to participants.  
Missing or extra copies are available from Room 229.


Q-SWITCH Nd:YAG LASERS  
IN OPHTHALMIC MICROSURGERY

- VAST POPULARITY ATTAINED BY THE Nd:YAG LASER IN OPHTHALMOLOGY FOR THE DISRUPTION OF TRANSPARENT MEMBRANES BOTH IN THE ANTERIOR AND IN THE POSTERIOR CHAMBER OF THE EYE
- Q-SWITCHING Vs MODE-LOCKING
- THIRD GENERATION OF Nd-YAG Q-SWITCH LASERS: COMPACTNESS AND RELIABILITY


DEVELOPMENT OF A Q-SWITCH  
Nd:YAG LASER FOR OPHTHALMIC  
APPLICATIONS AT C.E.Q.S.E-EVLASER

- THE REFLECTING CLOSE-COUPLED CAVITY: A COMBINATION OF HIGH EFFICIENCY AND EXCELLENT BEAM QUALITY
- STABLE AND UNSTABLE RESONATOR CONFIGURATIONS
- EFFICIENCY ENHANCEMENT WITH PLASTIC ENERGY CONVERTERS
- PLASTIC OR GLASS PASSIVE Q-SWITCH ELEMENTS FOR RELIABLE OPERATION


# COMPARISON OF DIFFERENT LASER PUMP CAVITIES:




# COMPARISON OF DIFFERENT LASER PUMP CAVITIES: BEAM DISTRIBUTION




# Nd:YAG LASER CONFIGURATIONS






a. unstable resonator



b. stable resonator



# INPUT-OUTPUT CURVES OF THE Nd:YAG LASER



## PERFORMANCES OF THE Nd:YAG LASER

- PUMP CAVITY: Reflecting close-coupled
- PUMP ENERGY: 0-30 J
- PUMPING EFFICIENCY: 2%
- OUTPUT ENERGY (Free-running)
 

ST	0-500	0-30
UN	0-400	0-80
- OUTPUT ENERGY (Q-Switch)
 


ST	300	15
UN	200	80
- DIVERGENCE
 

ST	4	1.2
UN	1.5	0.8
- STABILITY ( $\Delta E/E$ )
 

ST	3%	2%
UN	1.1%	1.1%
M.M.		S.M.

INPUT ENERGY (J)

# LAYOUT OF THE EXPERIMENTAL APPARATUS



## DETERMINATION OF BREAKDOWN THRESHOLDS IN LIQUIDS

- Laser pulse Energy: 0-80 mJ M.M.
- Laser pulse duration: 5 or 12 ns.
- Cone angle:  $11^\circ$  or  $16^\circ$   
(without B.E. for M.M., with B.E. for S.M.)

## WITH MULTIMODE Q-SWITCH LASER

