

SMR 1495 - 3

WINTER COLLEGE ON BIOPHOTONICS:
Optical Imaging and Manipulation of Molecules and Cells
(10 - 21 February 2003)


Fluorescence Lifetime Imaging

P.M.W. FRENCH
Head of Photonics Group, Physics Department,
Imperial College London, U.K.

These are preliminary lecture notes, intended only for distribution to participants.


Imaging tissue with fluorescence

Aim: to detect or image different types of tissue or states of tissue using optical radiation to achieve **contrast**



Problems: heterogeneity, scattering and background fluorescence

Difficult to make absolute intensity measurements




Solution: image $I(\lambda)$ or $I(\tau)$ – **relative measurements**

Laser-induced fluorescence (LIF) imaging of biological tissue

Fluorescence contrast

Most biological samples exhibit "autofluorescence"


Challenge is to find an unambiguous contrast between target and background, e.g. *normal and malignant tissue*



Tissue differences may be enhanced using a contrast agent, e.g. *HPD, ALA-S-induced Porphyrin IX*



Wavelength-ratiometric laser-induced fluorescence imaging of brain tissue



Wavelength ratiometric fluorescence measurements reduce the impact of scattering, since this is assumed to be constant with wavelength – or can be fitted to an assumed model

Experimental data from Lund, Svanberg et al.

Spectrally-resolved fluorescence imaging

How to acquire spectral data cube (x, y, λ) ?

Scanning microscopy $\Rightarrow$ serial pixel acquisition (single channel detection)
 $\Rightarrow$ slow image frame rate
 $\Rightarrow$ can use sophisticated detectors


e.g. spectrograph, spectrometer, FTS with PMT, APD detectors

Wide-field (x,y) imaging $\Rightarrow$ parallel pixel acquisition (many channels)
 $\Rightarrow$ high speed image acquisition
 $\Rightarrow \lambda$ requires one dimension $\Rightarrow$ need to scan?

- acquire x, λ , scan y : e.g. spectrograph with CCD detectors
- acquire x, y , scan λ : e.g. filter wheel, electro-optic filter with CCD detectors
- acquire x, y , scan z : e.g. FTS with CCD detectors

Wide-field multispectral fluorescence imaging

Wide-field imaging $\Rightarrow$ parallel pixel acquisition
 $\Rightarrow$ high speed data acquisition




Multi-spectral imager provides 2-channel (up to 8) spectral resolution
- preserves whole-field 2D/3D imaging

(Dichroic beamsplitter may be replaced by a polarizing beamsplitter to provide whole-field images of time-resolved polarization anisotropy)


Imaging tissue with fluorescence

Aim: to detect or image different types of tissue or states of tissue using optical radiation to achieve **contrast**



Problems: heterogeneity, scattering and background fluorescence


Difficult to make absolute intensity measurements



Solution: image $I(\lambda)$ or $I(\tau)$ – **relative measurements**

How to measure fluorescence lifetime?

Time-domain measurement




Frequency-domain measurement

$$\tan \phi = \omega \tau_p$$


$$m = [1 + \omega^2 \tau_m^2]^{-1/2}$$

For single exponential decay: $\tau_p = \tau_m = \tau$




Fluorescence lifetime imaging technology


Time domain



Frequency domain



Fluorescence lifetime imaging (FLIM)



Ultrafast laser sources for FLIM

Ti:sapphire femtosecond laser (Spectra-Physics)

100 fs @ 415 nm, 100 pJ @ 80 MHz

Diode-pumped ultrafast Cr:LiSAF laser oscillator and amplifier system (FLOAT)

10 ps @ 415 nm, 1 pJ @ single shot, 1 Hz-20 kHz

Blue diode laser (PicoQuant)

FWHM < 30 ps @ 400 nm,
Average Power = 1.1 mW,
Repetition rate up to 40 MHz



FLIM of dye samples: chemically specific imaging

Dye samples:

Coumarin 314	DASPI	Coumarin 314	DASPI	Coumarin 314
--------------	-------	--------------	-------	--------------

Fluorescence lifetime map:










FLIM of fluorophore environment (viscosity)

DASPI solvent: $\frac{\text{ethanol}}{\text{glycerol}}$


20	40	60	80
80	60	40	20

Fluorescence lifetime map:






Stretched Exponential Function



- 1863: F. Kohlrausch describes structural relaxation in glassy fibre.
- 1970: G. Williams & D.C. Watts describe dielectric relaxation in polymers.
- since then: applied to earthquakes, galactic light emissions, economics, biological extinction, MRI, etc.


Stretched Exponential Function



Plotting parameters


- average lifetime $\langle \tau \rangle = h \cdot \tau_{0.01} \cdot \Gamma\left(\frac{1}{h}\right)$
- heterogeneity $h \sim \text{width}$

Complex fluorescence decay profiles in rat ear (x10)



Laplace transformation: $I(t) \Rightarrow \tau, \rho(\tau)$

CONTIN algorithm by S.W. Provencher: *Comput. Phys. Commun.* 27, 213 (1982).



How to represent lifetime, heterogeneity and intensity information


H,S,L colour plot: H,S encode lifetime, L encodes intensity















H,S,L colour plot: H,S encode heterogeneity, L encodes intensity



Spectrally-resolved FLIM of rat heart tissue








Fluorescence Lifetime and Anisotropy Imaging (FLAIM): application to cells?

CFSE stained B cells



Outlook for microscopy and biomedical imaging

- More functional imaging
 - ... FLIM and FRET, multi-spectral imaging
 - ... Super-resolution
 - ... PSF engineering, e.g. 4 Pi microscope, STED
- Multi-modal microscopes
 - ... Combining confocal microscopy with OCT, THG, multi-photon fluorescence microscopy
- High-speed imaging
 - ... Wide-field optical sectioning
 - ... Multi-foci microscopes
 - ... High speed cameras and multi-channel imaging
- Compact user-friendly low-cost laser technology

For more information...

For a review of the development of ultrafast laser technology:

French, P. M. W., *The Generation of Ultrashort Laser Pulses: Reports on Progress in Physics*, 1995, 58(2): p. 169-252.

For a review of biomedical optics:

French, P. M. W., *Biomedical Optics in the 21st Century*, *Physica Scripta*, (June 1999) 41-48
<http://physicsweb.org/article/world/12/6/8>

For a review of fluorescence lifetime imaging:

Cole, M. J., J. Siegel, S. E. D. Webb, R. Jones, K. Downing, M. J. Dayel, D. Parsons-Karavassilis, P. M. W. French, M. J. Lever, L. O. D. Stedman, M. A. A. Neil, R. Juskaite, and T. Wilson, *Time-domain whole-field fluorescence lifetime imaging with optical sectioning*, *Journal of Microscopy*, 203 (2001) 248-257
 See also *OPN* November 2002 at www.csa.org

<http://photonics.ic.ac.uk>

Look for biomedical optics pages and recent conference presentations available online

Further reading

<http://micro.magnet.fsu.edu/primer/index.html>

Denk, W., Stuckler, J. H. & Webb, W. W. 2-Photon Laser Scanning Fluorescence Microscopy. *Science* 248, 73-76 (1990).

Müller, M., Squier, J., Wilson, K. R. & Brakenhoff, G. J. 3D microscopy of transparent objects using third-harmonic generation. *Journal of Microscopy-Oxford* 191, 266-274 (1978).

Webb, R. H. Confocal optical microscopy. *Reports on Progress in Physics* 59, 427-471 (1996).

Anderson Engels, S., Khairaberg, C., Svanberg, K. & Svanberg, S. In vivo fluorescence imaging for tissue diagnostics. *Physics in Medicine and Biology* 42, 815-824 (1997)

Cole, M. J. et al. Time-domain whole-field fluorescence lifetime imaging with optical sectioning. *Journal of Microscopy-Oxford* 203, 246-257 (2001)

Schrader, M. & Hell, S. W. 4Pi-confocal images with axial superresolution. *Journal of Microscopy-Oxford* 183, 189-195 (1996)

Straub, M. & Hell, S. W. Multifocal multiphoton microscopy: a fast and efficient tool for 3-D fluorescence imaging. *Bioimaging* 6, 177-185 (1998)

Klar, T. A., Jakobs, S., Dyba, M., Egnier, A. & Hell, S. W. Fluorescence microscopy with diffraction resolution barrier broken by stimulated emission. *Proceedings of the National Academy of Sciences of the United States of America* 97, 8206-8210 (2000)