

Comparative study of label-free electrochemical immunoassay on

various gold nanostructured electrode

S. Rafique and A. S. Bhatti Centre for Micro & Nano Devices, Department of Physics , Park Road,

COMSATS Institute of Information Technology, Islamabad, PAKISTAN.




- >In labeled immunosensor a 'label' is also incorporated into the assay to measurability.
- Most popular labels in immunoassays is Enzyme-linked immunosorbent assay (ELISA).
- >ELISAs includes multiple steps of incubation and washing.
- ≻It being an enzymatic reaction, even small
- quantities of non-specific binding might result in false signal

Advantages of label free detection


- > It eliminates the need for tags, dyes or specialized reagents required for assav development.
- > Simplifying assay design, minimizing liabilities created by the use of labels

Schematics of immunosensor preparation



Nanstructure s	HAuClO ₄ (mM)	HCLO₄ (M)	Voltage (V)	Time (min)
Pyramid	40	0.1	-0.08	2
Spherical	40	0.1	-0.2	2
Rod	4	0.1	-0.08	2

Optimization of experimental conditions




- The maximum current obtained tpH = 7.0.


 The 2 hr incubation time was selected.
 To ensure the prevention of saturation and for better performance of immunosensor.

Surface Morphology










Where as nano rods had an average diameter about 120 nm.

Nanostructure s	Surface area (µm²)	Surface Coverage (%)
Pyramid	45.5	52.9
Spherical	42.3	63.9
Rod	52.1	40.7

FESEM and AFM images of (a) Pyramid (b) spherical (c) Rod like

Cyclic Voltammetry



- > Peak current of the bare Au was 8.69µA.

 The peak current increased
 1.40 times as compared to bare Au
- electrode.

 > With further modification with
- anti CEA showed 0.83, 0.76 and 0.77
 times drop in the peak current for
 pyramid, spherical and rod like nanostructures, respectively.

Limit of Detection (LOD)




- The current dropped with the increase in the concentration from 1pg/ml to 1000ng/ml.

 The limit of detection (LOD) was calculated using the

Nanostructures	Bare Au	Pyramid	Spherical	Rod
LOD (ng/ml)	0.07	0.0039	0.0036	0.0045
Sensitivity (µA ng ⁻¹ . ml)	0.211	0.339	0.457	0.312


Association constant



The values came out to be 0.0695, 0.0707, 0.0948 and 0.0696 for bare Au, pyramid, spherical and rod like nanostructures , respectively.

Specificity and stability of immunosensor





Conclusion

- The spherical nanostructures were smaller in size and has larger value of surface area as compared to the
- pyramid and rod-like nanostructures.

 Due to the higher surface area the spherical nanostructure showed better electrochemical performance than the other types of nanostructures.
 - It was found the limit of detection of the nanostructured electrodes was 4 pg/ml and have stability almost for a month

Reference

- T. P. Sai A. K. Raychaudri, J. Phys. D: Appl. Phys. 40 (2007) 3182.
- G. legay, L. Markey et al, Ultramicroscopy, 107 (2007) 1111.
- G. Demirel, B. Garipean et al, Surface Science, 602 (2008) 952.
- F. W. Delrioo, C. Jaye et al, Applied Physics Letters 94, (2009) 131909. S. Sampath, Indiand Society of Sci 65 (2005) 753.

Acknowledgement

- Higher Education Commission (HEC), Pakistan for
- All ASB Lab fellows